The science behind Firstbeat HRV analytics

Tero Myllymäki Head of physiology



AUTONOMIC NERVOUS SYSTEM CONTROLS OUR BODY

Parasympathetic





Sympathetic



Adapted from: OpenStax College - Anatomy & Physiology, Connexions Web site



HRV GIVES US ACCESS TO AUTONOMIC NERVOUS SYSTEM

Heart rate variability (HRV) refers to beat-by-beat changes in heart rate.

HRV is mediated by **autonomic nervous system** (ANS) ¹ and scientific studies have shown HRV can be used as a non-invasive measure of sympathetic and parasympathetic activity.

HRV reflects the control of heart. Heart reflects body's physiological processes and demands such as exercise status, hormonal reactions, metabolism, cognitive processes, stress reactions, relaxation/recovery, sleeping and emotions.



Autonomic nervous system (ANS): Balance, Control of Heartbeat

+ Sympathetic Activation, mobilization of body resources

Parasympathetic Calming down, recovery, restoration of resources

¹ Task Force 1996. Heart rate variability: standards of measurement, physiological interpretation and ^{20 May 2019} clinical use. Circulation. 93(5):1043-65.

WHERE IS THE PROOF?

HR & HRV after ANS blockade



Ability of short-time Fourier transform method to detect transient changes in vagal effects on hearts: a pharmacological blocking study. Am J Physiol Heart Circ Physiol 290: H2582–H2589, 2006. Kaisu Martinmäki, Heikki Rusko, Sami Saalasti, and Joni Kettunen.



HRV AND RESPIRATION

RR-intervals (ms)



Time



HEALTHY HEART HAS VARIATION



Reduced / low HRV is associated with

Acute stress (Hall et al 2004)

Work stress (Vrijkotte et al 2000, Clays et al 2011)

Heart disease, anxiety, depression, asthma and PTSD ...

High HRV is associated with

Reduced morbidity and mortality (Sajadieh et al 2004; Stein et al 2005) Psychological well-being and quality of life (Geisler et al 2010) Better physical fitness (de Meersman 1993)



CONCEPTUAL FRAMEWORK OF VAGAL CONTROL AND HRV

Unifying conceptual framework of factors influencing cardiac vagal control



Laborde S, Mosley E & , Mertgen A (2018). A unifying conceptual framework of factors associated to cardiac vagal control. Heliyon 4 (2018) e01002. doi: 10.1016/j.heliyon.e01002



THE DIFFICULTY





OUR APPROACH

Firstbeat wants to truly understand how human body works for **providing the most accurate and meaningful feedback** for anyone interested on wellbeing, health and performance.

Firstbeat analytics is **based on applying heart rate variability** and other sensor information and packing that into compelling UX on human physiology.

Firstbeat analytics has been developed to form a comprehensive network of algorithms during 20 years **by expert teams** of physiology, mathematics, and software development.





Recovered

7h 40min 22:45 – 6:25

Calm sleep

1 2

Sleep period was long enough and calm. You're well recovered.

ANALYTICS IS A MUST FOR INSIGHTS





DIGITALIZING BODY REACTIONS

- 1. Measurement (RR-intervals, acceleration, speed, power)
- 2. Signal processing, artefact correction, and quality control
- 3. Algorithms, decision trees, neural networks, Al
- 4. Digital modelling of physiology systems
- 5. Personal calibration and adaptation

20 May 2019

- 6. Body status, interpretation for different use cases, insights on the contexts
- Benefiting from the knowledge, e.g. feedback on lifestyles and effects of behaviors

Image from: "Anatomy in Motion." Authored by: Beth Scupham. Located at: <u>https://flic.kr/p/c98TPd</u> Content Type: CC Licensed Content, Shared Previously. License: CC BY: Attribution.





FROM LABORATORY TO REAL-LIFE

- User profile
- Heart rate variability
- Heart rate level
- Electrocardiography (ECG). Tidal volume
- Optical heart rate
- Optical pulse waves
- Acceleration
- Step rate / cadence
- GPS or running speed
- Cycling watts
- Altitude
- Temperature
- Duration
- Distance
- Oxygen uptake
- Excess Post-Exercise Oxygen Consumption (EPOC)
- Energy expenditure

- VO2Max
- Respiratory rate
- Ventilation
- Percentage of vital capacity
- VO2 and VCO2 concentration
- Respiratory exchange ratio
- Blood lactate
- Neuromuscular parameters
- Cortisol and other blood markers
- Electroencephalography (EEG)
- Sleep stages
- Perceived exertion
- Perceived stress
- Questionnaire data
- Activities / real-life simulations
- Journal / diary / training logging

FIRSTBEAT

FIRSTBEAT ANALYSIS DEVELOPMENT: VO2 AND VO2MAX



HR



FIRSTBEAT ANALYSIS DEVELOPMENT: VO2MAX EVOLUTION



APPROACHES ON DEVELOPMENT AND VALIDATION OF STRESS

Biogenic: ANS blocking

• <u>CAUSE OR BLOCK stress</u> <u>response</u> directly.



Psychosocial: stress test

 <u>Manipulate condition</u>, subject people to stressful or relaxing condition



Real-life setting

 <u>Observe</u> people during reallife, measure physiological stress & subjective stress



High to very high

Moderate to high

Moderate to low

Expected max correspondence with HRV measured stress

20 May 2019



TRIER SOCIAL STRESS TEST (TSST) STRESS TEST



A.P. Allen et al. (2014). Biological and psychological markers of stress in humans: Focus on the Trier Social Stress Test. Neuroscience and Biobehavioral Reviews 38: 94–124



20 May 2019

EXAMPLE OF STRESS TEST DATA



- Baseline rest 10 min
- Speech preparation 10min
- **Public speech** 5 min (in front of 2 judges, videotaped)
- Mental arithmetic 5 min (counting backwards from number 1022 in steps of 13)
- **Rest** 10 min
- HRV measurements
- Cortisol measurements
- Perceived stress



EXAMPLE OF FIELD STUDIES ON STRESS

The stress generated or maintained outside working hours correlates significantly with a lower quality of recovery during the 24 h workday. It is necessary to prioritize strategies that help improve stress management in executives through the improvement of tools and strategies that promote greater relaxation outside working hours mainly.



Crespo-Ruiz et al. 2018, Executive Stress Management: Physiological Load of Stress and Recovery in Executives on Workdays. Int. J. Environ. Res. Public Health, 15, 2847.

Environmental Research and Public Health	MDPI
Article Executive Stress Management: Pl Stress and Recovery in Executive	hysiological Load of s on Workdays
Beatriz Crespo-Ruiz ^{1,2,*} , Shai Rivas-Galan ² , Cristina Fern and Luis Maicas-Perez ¹	andez-Vega ² , Carmen Crespo-Ruiz ²
¹ Physical Activity and Sport Education Department, Faculty of Castilla-La Mancha, 45002 Toledo, Spain; info@fluismaicas.com 2 Freedom & Elow SL, 28000 Madrid, Spain; shair/was@freedom Cristina.fernandez@freedomandflowcomp.any.com (CC-RV); carmer.newspo@freedomandflowcomp.any.com (CC-R) - 1 Correspondence: beatriz.crespo@fuclmes; fic1: +34-617-41511	Sport Sciences of Toledo at University of andflowcompany.com (S.RG.);
Received: 8 October 2018; Accepted: 2 December 2018; Published:	13 December 2018 check for updates
Abstract: Objective: The use of high-performance sports techn of stress and the quality of recovery in a population of exect Heart rate variability values were recorded during 48 h stress/recovery quality (stress balance) was obtained for th work, and night in a workday. <i>Results:</i> We observed a neg measurement in the course of a workday, being negative e night. The stress generated or maintained outside working h quality of necovery during the 24 h workday. <i>Conclusions:</i> It help improve stress management in executives through the mainly promote greater relaxation outside working hours. Keywords: stress; stress management; human resources; ex- business	nology to describe the physiological load trives during the workday. Methodology: from which the relationship between ee differentiated time slots: work, after ative stress balance during the 24 h of t work and after work, and positive at ours correlates significantly with a lower is necessary to prioritize strategies that mprovement of tools and strategies that ecutives; physiology; health; technology;
1. Introduction	
Stress is defined as the physical and mental responses or perceived changes and challenges in life. Stress is, for many 21st century. However, it was not officially classified as a dis differently to anxiety until the 5th edition of the Diagnostic an (DSM) by the American Psychiatric Association published or Newadays society faces the rerest challenge of adapting to	of the body and adaptations to real and , considered as the silent disease of the sorder or alteration to take into account d Statistical Manual of Mental Disorders 1 October 2016.
to ormany society nates use great channenge of adapting it to a large number of internal and external stimuli that, poord stress triggers, with the physical and mental consequences working time. To this effect, there are several challenges that when it comes to promoting health and well-being. Actually it arxiely, and/or depression, have led to financial losses up percentages of more than 75% of workers in the European Un more than 20% the working age population that suffer serio stress syndrome known as burnout [1]. In this sense, and due to the multifactorial nature of stress	a commodus rate of change, connected managed, can become important daily that all of this involves in and out of companies and governments must face the increase of pathologies, such as stress, to €136 million in companies reaching ion who suffer daily stress being already us health problems, such as the chronic , other factors such as poor rest, physical
inactivity, overweight, and obesity are increasingly present	in today's companies, aggravating and



FIRSTBEAT ANALYSIS DEVELOPMENT: SLEEP



MMMMMM

NREM slow wave sleep (SWS) REM sleep



CONQUERING THE HRV SUMMIT. WHERE IS IT?





LAYERS OF SCIENCE APPROACH





Thank you!

Tero Myllymäki tero.myllymaki@firstbeat.com

